

ANSI/AAMI ST91: 2015

Flexible Endoscope Reprocessing and the Importance of AAMI ST91




April 21, 2017

Fred Alston, CSPDT
Director of Clinical Sales and Services
Healthmark Industries

Objectives

- Discuss the key provisions and competency recommendations of the standard
- To identify best practices in reprocessing of flexible endoscopes
- Discuss methods of cleaning verification and surveillance testing to determine if an endoscope is patient ready

What is ANSI/AAMI ST 91?



- Flexible and semi-rigid endoscope reprocessing in health care facilities
- Contains best practices for scope reprocessing in ANY setting
- Available for purchase at the www.aami.org

Risk of Endoscopy Related infection or Other Adverse Patient Reactions

- Spread on infections related to endoscopy:
 - **Exogenous** infections = Microorganisms spread from patient to patient by contaminated or malfunctioning scopes or equipment
 - Microorganisms may be transmitted from patients to endoscopy personnel and/or from endoscopy personnel to patients
 - **Endogenous** infections = Microorganisms spread from the GI tract through the bloodstream during an endoscopy procedure to susceptible organs, or may spread to adjacent tissues that are breached as a result of the endoscopic procedure

Risk of Endoscopy Related infection or Other Adverse Patient Reactions

- Other risks related to endoscopy:
 - Chemical substances can remain on devices from various chemicals used during the procedure or processing that can cause toxic reactions in subsequent patients.
 - Chemical burns, colitis, anaphylaxis, death
 - Devices may be damaged or rendered difficult to use due to mishandling or inadequate processing.

Objective – ST91

- Provide guidelines for processing of flexible endoscopes
 - Includes all stages of reprocessing HLD and sterilization of scopes and accessories
- Include flexible gastrointestinal (GI) endoscopes; bronchoscopes; ENT scopes; surgical flexible endoscopes (e.g., ureteroscopes); and semi-rigid operative scopes (e.g., choledochoscopes)
- Exclusions
 - Rigid endoscopes and probes (e.g., TEE probes)



ST91 Scope – What's contained in this standard?

- Definitions
- Design of endoscope processing areas
- Personnel considerations
- Cleaning
- High level disinfection
- Automated endoscope reprocessors (AERS)
- Liquid chemical sterilization
- Gaseous chemical sterilization
- Processing accessories
- Storage and Transportation to site of use
- Quality Control including cleaning verification
- Quality Process Improvement
- Informational Annexes

Best practices for processing flexible endoscopes

- Meticulous attention to all steps in processing endoscopes, their components and accessories is critical making them safe for subsequent patient use
- Steps are outlined in the document in detail and include the following categories
 - Precleaning, transportation, leak testing, cleaning, rinsing, inspection or testing for cleanliness, high-level disinfection & sterilization and monitoring of the process, rinsing, drying, alcohol flush, & storage

Highlights of AAMI ST 91

- Gives recommendations for:
 - **Certifications for technicians performing reprocessing**
 - **Monitoring the manual cleaning process**
 - **Monitoring the automatic cleaning process**
 - **Monitor water quality**
 - **Monitor temperature**
 - **After cleaning, all detachable valves should be kept together with the same endoscope as a unique set**
- Risk Assessment
- Proper documentation and quality assurance parameters




Processing / Reprocessing

Processing (or reprocessing) is a process carried out on a device to allow its subsequent safe use, which can include cleaning, disinfection, sterilization, and related procedures



Best practices in Precleaning

- Prevents buildup of bioburden, development of biofilms, drying of patient secretions
- Occurs at point of use immediately after the procedure
- Don fresh PPE
- Prepare a cleaning solution (or water if validated) according to the solution manufacturer's written IFU.
- Wipe insertion tube with a low or non-linting cloth/sponge soaked in the freshly prepared cleaning solution.
 - Note: cloth/sponge is single-use only



Remember to follow the IFU for the endoscope and detergent!

Best practices in precleaning

- Ensure that controls are in the free/unlocked position.
- Suction solution through the suction channel as per manufacturer's written IFU.
- Flush the air/water channels with solution using the cleaning adapter per manufacturer's IFU.
- Flush all other channels (e.g., auxiliary water or elevator channels) with solution, if present.
- Suction the solution through the endoscope until clear.
- Detach the endoscope from the light source and suction pump.
- If applicable, attach the fluid-resistant cap.
- Visually inspect the endoscope for damage.

Contaminated Transport

- From procedure room to reprocessing area:
 - Closed, labeled transport containers
- Place a single endoscope in a container by naturally coiling it in large loops.
- Separate endoscopy accessories from the endoscope to prevent puncture and damage.
- Labelled appropriately as biohazard



Best practices for Leak Testing

- Occurs in processing area prior to immersion in cleaning solution.
- Serves to detect damage that would allow for fluid-invasion
- Wear PPE
- Ensure fluid-resistant cap is on prior to submersion
- Use a basin of **water** or surface large enough to ensure that the endoscope is not coiled too tightly to mask holes.
- Allow for sufficient time to observe the endoscope for leaks, manipulate knobs and buttons



Best practices for Leak Testing

- Outlines 4 general methods for performing leak test:
 - Manual (dry) leak testing
 - Mechanical (wet) leak testing
 - Mechanical (dry) leak testing
 - Mechanical AER leak testing
- Refer to manufacturer's IFU for detailed steps
- For failures, refer to manufacturer's IFU for modified processing steps being sure to maintain positive pressure throughout



Best practices for manual cleaning

- Soil remaining on the endoscope may interfere with the ability of the disinfection or sterilization process to effectively kill or inactivate microorganisms
- If process is not initiated immediately, follow written IFU for delayed reprocessing from manufacturer
- General process is outlined including
 - Don fresh PPE, use fresh detergent solution, monitor the temperature of the cleaning solution


Best practices for manual cleaning

- Cleaning steps:
 - Clean with a single-use lint-free cloth/sponge
 - Submerge scope to prevent splashing contaminated fluids
 - Use a cleaning brush with specifications per manufacturer's IFU
 - Brush all channels, cylinders, openings and forceps elevators per IFU



Best practices for manual cleaning

- Cleaning steps (continued):
 - Use recommended cleaning adapters
 - Flush all channels, rinse all channels, air purge all channels
 - Repeat until there is no visible debris
 - Soak, scrub, brush & rinse all reusable/removable parts
 - Automated flushing pumps may be used during manual cleaning



Cleaning Solutions (Detergents)

- Designed for endoscope cleaning
- Typically neutral detergents
 - May or may not contain enzymes
 - Numerous products available
- Essential features
 - Optimum cleaning performance
 - Manufacturers' labeling
 - Device protection
 - Water quality control
 - Toxicity validation



Automated flushing systems

- If a flushing pump is used, follow manufacturer's written IFU
- Ensure compatibility of endoscope with model of flushing system
- Use fresh solution with each endoscope
- Clean and disinfect tubing and equipment according to manufacturer's IFU
- Perform any other QA testing as recommended (e.g. daily volume verification)



Rinsing after cleaning

- Thoroughly rinse with copious volumes of potable water
 - AAMI TIR34
- Follow IFU of endoscope & cleaning solution to determine the amount of water needed for rinsing, psi/pressure, and number of rinses
- Use recommended cleaning adapters
- Rinse all external and internal surfaces
- Perform an air purge of all channels
- Dry exterior with a lint-free cloth/sponge
- **Keep detachable valves together with the same endoscope as a unique set**



Best Practices for Cleaning Verification & Process Monitoring

- Cleaning verification is performed following cleaning to verify the effectiveness of a cleaning process **PRIOR TO DISINFECTION**
- Cleaning verification should include:
 - Visual inspection
 - Testing of the cleaning efficacy of mechanical equipment
 - Monitoring of key cleaning parameters
- Use of methods to detect organic residue should be considered

Verifying Clean through Inspection

- Visual inspections and testing of the equipment
 - Inspecting organic residues
 - Testing for any cracks in the devices
 - Checking integrity of fiber optic bundles
- Use lighted magnification and inspect throughout process
- Consider inspection with borescope
 - ST91 and AORN recommendations
- SGNA – Treat as a **“Time out”**
- Methods to measure organic & other residues found on scopes
 - Protein
 - Hemoglobin
 - Carbohydrates
 - ATP



Optical & Enhanced Inspection

AORN Recommendations:


- Visually inspect with lighted magnification for cleanliness, integrity, and function before use, during the procedure, after the procedure, after cleaning, and before disinfection or sterilization.
- Inspection helps to identify residual organic material and defective items and remove from service soiled/defective items that might put patients at risk for infection or injury.
- An endoscope that appears clean may harbor debris that cannot be seen without magnification. Lighted magnification may increase the ability to identify residual soil or damage.
- Internal channels of endoscopes may be inspected using a borescope. Borescopes penetrate the lumen and allow for improved visual inspection.



Cleaning verification recommendations

- Current recommendations support testing of the manual cleaning process at pre-established regular intervals:
 - AAMI ST91: Regular intervals, i.e. **Weekly or preferably daily**
 - AORN: Regular intervals such as with **EACH reprocessing cycle** or daily
 - SGNA: Confirm the adequacy of manual cleaning by using a rapid cleaning monitor. If the tool results are positive, this allows for the re-cleaning of the endoscope prior to disinfection. **Frequency determined by facility.**

Manual Cleaning Verification Monitors



Which Organic Parameters to monitor?

Flexible endoscope biopsy channel: (Alfa et al 2002)

- Protein; $< 6.4 \mu\text{g}/\text{cm}^2$
- Carbohydrate; $< 1.8 \mu\text{g}/\text{cm}^2$
- Hemoglobin; $< 2.2 \mu\text{g}/\text{cm}^2$
- Endotoxin; $< 2.2 \text{ EU}/\text{cm}^2$




Best practices for High-Level Disinfection

- Standard of care for reprocessing semi-critical instruments
 - Those devices which contact mucous membranes
 - Sterilization preferred or HLD with an FDA-cleared HLD prior to next use
- HLD defined as a germicide that inactivates all microbial pathogens, except large numbers of bacteria endospores when use in accordance with labeling

Best practices for High-Level Disinfection

- HLD types include:
 - Glutaraldehyde
 - OPA
 - Peracetic Acid
 - Chlorine
 - Hydrogen Peroxide
 - Combination products
- Refer to www.fda.gov for a list of cleared HLDs
- HLD can be performed manually or with an automated endoscope reprocessor (AER)

High-level Disinfection



Best practices for High-Level Disinfection


- Reusable HLDs must be monitored to ensure that it is above the Minimum recommended concentration (MRC)
 - Test prior to each use per IFU
 - Solution is used repeatedly until it fails test strip or meets its maximum use life, whichever comes first
 - Do not “top off” HLD unless instructed by HLD manufacturer
 - Can not be used to extend the use life of HLD

Best practices for High-Level Disinfection

- Single-use HLD's:
 - Used with specific AER's
 - Can be concentrated or ready to use
 - Examples concentrated OPA and peracetic acid
 - MRC is tested through either test strip or chemical monitoring by the AER
- HLD's need to contact ALL surfaces
 - Internal channels and external surfaces
 - Complete immersion
 - Monitor exposure times precisely
 - Remove air bubbles from surfaces of endoscope

Remember to Rinse!

- Rinsing is often overlooked and underestimated
 - Removal of chemicals and residual soil such as protein (e.g., enzymes used during cleaning)
 - Devices should not present a toxic risk to patients
 - Water quality/purity can impact this
 - Number of rinses and rinsing method using fresh water with each rinse



Key Points with Disinfection

- Label claims can vary
 - Safety, preparation, contact time, numbers or rinsing etc
 - Request specifics from manufacturers (e.g., 'rinse thoroughly')
- Single use or multiple use
 - Use of solution test strips to verify minimum recommended concentration (MRC)
- Multiple use disinfectants
 - Closely reuse label claims, including maximum reuse life
 - 'Topping off' of solutions
- All parts of the device should be contacted
- Importance of rinsing
 - Correct water quality (bacteria-free; AAMI TIR34)
 - Fresh water for every rinse (by immersion)
 - Correct number of rinses
- Device inspection prior to use



Use of Automated Endoscope Reprocessors (AER)

- Machines designed to clean and/or disinfect endoscope and components using an LCS/HLD solution
- Use of AER's may be more efficient and leads to less user exposure and helps to ensure repeatable results
- Section has detail on types of AER available and features of their cycles
- If AER cycle is interrupted, it should be repeated
- Purchase considerations are outlined

Automated Endoscope Reprocessors (AERs)

- Review claims and instructions for use carefully
 - FDA clearance
 - Note any limitations such as disinfection of connector contact sites on endoscopes, device preparation, flushing capabilities for all lumens, drying capabilities etc
- Correct use of cleaning detergents (when applicable) and disinfectants
- Rinse water control
- Routine maintenance



Manual Drying and Alcohol Flush

- Effective drying reduces the risk of microbial contamination post HLD
- Waterborne organisms can pose an infection control risk to some patients
 - Bronchoscopy and ERCP patients
- Presence of such organisms in conjunction with retained moisture can lead to biofilms and patient risk
 - Especially true if tap water is used for final rinse
- Hanging to dry' or 'drip dry' is NOT effective
 - Most AERs 'purge' water from the endoscope lumens not 'dry'



Manual Drying and Alcohol Flush

- Drying is achieved by flowing air through the endoscope channels
 - Facilitate drying with alcohol flush (70-80% ethyl or isopropyl alcohol)
 - Follow endoscope IFU for amount to be used
 - **Follow with instrument quality forced air to ensure residual alcohol is removed**
 - Refer to endoscope IFU for psi recommendations
 - Dry all removable parts and do not reattach
 - Keep valves with the endoscope to ensure traceability




Liquid Chemical Sterilization (LCS)

- Liquid chemical sterilization system is used for heat-sensitive, critical medical devices when traditional methods are not feasible or available
- Devices are treated with LCS & rinsed with water
- Rinse water is treated but may not be not sterile
- Can not maintain sterility
 - Immediate use, no storage
- System is currently available
- Follow written IFU's of LCS system for proper use



Best practices for Sterilization

- Sterilization processes for flexible and semi-rigid scopes is discussed in detail
- Recommended for devices entering sterile body cavities
- Section outlines special considerations for terminal sterilization with primary source of info being endoscope's IFU
- More modalities compatible with surgical flexible endoscopes



Best practices for Sterilization

- Packaging considerations are outlined
 - Pouches, wraps, rigid containers
- Guidance for different sterilization modalities is given:
 - Steam
 - Ethylene Oxide
 - Hydrogen Peroxide gas
 - Ozone
- Sterilization is dependent on adequate **cleaning**, rinsing and device preparation
 - Drying is also essential
 - Packaging requirements (if applicable)



Storage of reprocessed endoscopes

- Endoscope should be hung vertically with the distal tip hanging freely in a well-ventilated, clean area following endoscope manufacturer's IFU for storage
- Angulation locks in the free position
- Sufficient space between endoscopes
- All removable parts should be detached, but kept together with the endoscope
 - (small bag or similar device)



Storage of reprocessed endoscopes

- Have policies and procedures in place regarding storage
- Ensure that endoscopes are adequately dry prior to placing in storage to prevent bacterial growth and biofilm
- Sterilized endoscopes should be stored in their container or packing in which they were sterilized



TRANSPORTATION IDENTIFICATION TAG		TRANSPORTATION IDENTIFICATION TAG	
DATE		DATE	
CASE PROCEEDURE		CASE PROCEEDURE	
ROOM/FACILITY		ROOM/FACILITY	
INITIALS		INITIALS	
CLEAN	PULL HERE	DIRTY	PULL HERE
TRANSPORT CHECKLIST		TRANSPORT CHECKLIST	
ALL INSTRUMENTS DISINFECTED	YES NO	ALL INSTRUMENTS DISINFECTED	YES NO
ALL INSTRUMENTS DRY	YES NO	ALL INSTRUMENTS DRY	YES NO
ALL INSTRUMENTS PACKED	YES NO	ALL INSTRUMENTS PACKED	YES NO
ALL INSTRUMENTS LABELED	YES NO	ALL INSTRUMENTS LABELED	YES NO
ALL INSTRUMENTS SECURED	YES NO	ALL INSTRUMENTS SECURED	YES NO
INITIALS		INITIALS	


Storage

- General considerations
 - Prevent coiling or kinking (hanging preferred)
 - Closed cabinets recommended
 - Tracking and traceability
- Hang time
 - Importance of risk assessment (facility-specific)
 - Policy and procedure development
- Liquid chemical disinfection or sterilization
 - Drying is essential
 - Reduce risks of recontamination
 - Transportation to point of use
- Gaseous sterilization
 - Correct storage conditions/methods



Risk assessment recommendations

- Risk assessment should be performed to address **length of storage (hang-time)**
 - Considerations should be given to the following:
 - Complexity of instrument, condition after processing (wet/dry, alcohol flush), transportation methods, conditions of storage environment, handling during storage, manufacturer's recommendations for storage, professional society guidelines, current research studies, protective devices to prevent
 - **Now in alignment with AORN recommendations to conduct a risk assessment**
- Develop protocols to ensure that users can readily identify an endoscope that has been processed and is ready for patient use.




Current recommendations for length of storage “hang time”

- AAMI ST91: Due to lack of consensus it is recommended to perform a **risk assessment** to establish maximum length of storage
- AORN: Perform a **risk assessment** with a multi-disciplinary team to establish a policy for maximum storage time that processed flexible endoscopes are considered safe to use without reprocessing.
- SGNA: 7 days based on a systematic review, if scopes are effectively reprocessed and stored in a way that keeps them completely dry and free from environmental and human contamination

Use of sterile endoscope sheaths

- Available for use with specified endoscopes
- Instructions for endoscope reprocessing for cleared devices recommend alternative processing instructions when the sheath remains intact after endoscope use
- Sheaths are not cleared for all types of endoscopes
- Two categories of sheaths:
 - Those intended to reduce the level of soiling of the endoscope
 - Those intended to prevent endoscope soiling and serve as microbial barriers
 - There are different processing instructions dependent on which type is being used
 - Refer to sheath manufacturer IFU and endoscope IFU



Processing of endoscope accessories

- Processing of certain endoscope components (valves) requires the same level of inspection, cleaning, and HLD or sterilization as the endoscope itself
- Process for manual cleaning & HLD of accessories is outlined
- Water bottle processing should occur according to manufacturer's IFU and at least daily

ST91

...

Education, Training and
Competency Recommendations

Education, Training and Competency

Recommendations

- All personnel performing processing be **certified** as a condition of employment.
- At a minimum, personnel should complete a certification exam.
- Frequencies of training/competency:
 - initial hire; annually; at designated intervals; or whenever new endoscopic models, new processing equipment, or products such as new chemicals are introduced for processing.

Education, Training and Competency

Recommendations

- Processing activities should be closely supervised until competency is verified for each processing task.
- Identification of items that are single-use
- Facility policies regarding sterilization and high-level disinfection, infection prevention, attire, hand hygiene, and compliance with local, state, and/or federal regulations

What Should The Education Include?

- Procedures for cleaning, leak testing, disinfecting or sterilizing, packaging, and storing each specific endoscope make and model, including equipment connections
- All aspects of decontamination (e.g., disassembly, manual and mechanical cleaning methods and how to monitor their effectiveness, microbiocidal processes, equipment operation, standard precautions, and engineering and work practice controls)
- Documentation of quality monitoring results.

What Should The Education Include?

- The operation of the specific cleaning processes & equipment, high-level disinfection processes, sterilizing systems and the methods used to verify operation
- Workplace safety, including OSHA standards for chemical use and biological hazards as well as workplace safety processes and procedures related to endoscope processing, high-level disinfection, and sterilization.

Quality Control Procedures

- Quality control is critical within endoscope reprocessing procedures
- Topics covered are product identification, traceability, documentation, record-keeping, verification and monitoring of HLD and sterilization process, product recalls and quality process improvement
- Facilities should develop comprehensive quality assurance and safety programs
- Each section outlines what parameters should be documented, tested, and/or maintained

Product Recalls

- Written policies should be in place for a recall event (HLD or sterilization failure)
- Policies developed in cooperation with infection prevention and risk management
- Establishing recall procedures helps to ensure patient safety, compliance with user facility reporting requirements to the FDA & allows for adequate follow-up actions




Microbial Surveillance

- Options include:
 - Traditional culturing
 - Gram negative test kits
- AAMI - No recommendation is made in the current version because of the timing of release.
 - Studies have identified the nature of microbial contamination likely to be found in improperly reprocessed endoscopes and have demonstrated the value of surveillance testing
- AORN: Base decision on a risk assessment
- Not ATP or cleaning verification tests



Guidance on culturing

- CDC Interim Guidance on culturing duodenoscopes updated 4/3/15
 - Sites to be cultured?
 - Instrument channel (suction/biopsy channel)
 - Distal end (elevators mechanism, elevator recess)
 - Elevator channel (on older, unsealed)
 - **Frequency: Every 30 days or 60 cycles**
- Mail back service for endoscope samples are now available



Monitoring for Gram-negative organisms in reprocessed scopes

- Enzymes specific to Gram-negative bacteria hydrolyze the substrate in the reagent vial
 - This generates fluorescence, which is read by the fluorometer, which then gives a reading.
- ST91: Types of verification testing may include enzyme based tests
 - Such as the gram negative test kit



Summary

- With heightened public concern and documented cases of improper reprocessing endoscopes, it is imperative that we must reducing the risk of exposure to improperly reprocessed medical devices.
- This is a shared responsibility among the healthcare facilities responsible for cleaning, disinfecting or sterilizing the devices.
- ST91 is your go-to guide for national standards in endoscope reprocessing and highlights best practices and quality control measures for each step along the way. Available at www.aami.org

My contact information:

Fred Alston

Cell: 724-448-4121

falston@hmark.com

www.hmark.com

www.hmarkgi.com